

# CREATININE "Mod. Jaffe"

Diagnostic reagent for determination of Creatinine concentration.

Liquid. Dual Reagents. Store at 15°C - 25°C

For in Vitro Diagnostic Use. Do not freeze

Ref No	Pack	Ref No	Pack	Ref No	Pack	Ref No	Pack
A2160	4X100 ml	T2160	4091 Tests	R2161	962 Tests	L2160	3000 Tests
A2161	4x50 ml	T2161	2841 Tests	BY2160	11136 Tests	L2161	1600 Tests
A2162	4x25 ml	S2160	2338 Tests	BY2161	6264 Tests	M2160	848 Tests
DM2160	960 Tests	S2161	1559 Tests	K2160	7000 Tests	M2161	1061 Tests
N2160	800 Tests	R2160	2885 Tests	K2161	5000 Tests		

## INTENDED USE

Test for quantitative determination of creatinin in human serum, plasma and urine.

## TEST PRINCIPLE

Creatinine reacts with picric acid in alkaline environment to form a color complex. Developing of this red color may be followed by photometrically at 500-520 nm. The association on surfactant and sodium tetraborate keeps interferences at minimum.

## TEST PARAMETERS

Method : Colorimetric, Kinetic, "mod."Jaffe, Increasing Reaction  
 Wavelength : Hg 510 nm (500 nm - 520 nm)  
 Temperature : 37°C  
 Sample : Serum, heparinized Plasma, Urine  
 Linearity : 0.20 mg/dL - 20 mg/dL

## REAGENT COMPOSITION

### Reagent 1: (R1)

Carbonate Buffer ≤ 120 mmol/L

Sodium Hydroxide ≤ 360 mmol/L

### Reagent 2: (R2)

Picric Acid ≤ 7.8 mmol/L

## REAGENT PREPARATION

Mix 1 unit of reagent 1 and 1 unit of reagent 2. Working reagents are stable up to 30 days at 15-25°C, well capped and away from light sources. Stability of unmixed reagents: up to expiration date on labels at 15-25°C; Stability since first opening of vials of unmixed reagents: ≥ 60 days at 15-25°C. For manual working procedures; if working reagent will be used; shake the Reagent 2 vial gently before pouring its contents into the Reagent 1 bottle. It is advisable to wash the Reagent 2 vial with a small volume of the prepared mixture in order to

completely rinse the vial and avoid any losses.

## REAGENT STABILITY AND STORAGE

Store at 2-8°C.

Reagents are stable until the expiry date shown on the label when stored tightly closed and if contaminations are prevented during their use.

Once opened vials are stable minimum 30 days at 2-8 C at optimum conditions. On board stability is strongly related to auto analysers cooling specification and carry-over values.

## SAMPLE

Serum, plasma, Urine. Creatinine is stable 24 hours at 2-8°C. Freeze samples for prolonged storage. Dilute urine sample 1:100 with deionized water. It could be convenient a slight acidification of urine with HCl.

## TEST PROCEDURE

### Sample Start

There have many ready application procedures dedicated to different kind of photometers and ready manual working process can be supplied on request.

There have many ready application procedures dedicated to different kind of biochemistry auto analyzers can be supplied on request.

### Substrate Start

There have many ready application procedures dedicated to different kind of biochemistry auto analyzers can be supplied on request.

## CALCULATION

$\frac{\Delta A \text{ Sample}}{\Delta A \text{ Standard}} \times \text{Concentration of Standard} =$   
 mg/dL (μmol/L) Creatinine in sample. For urine multiply result by 20.

**Unit Conversion**

mg/dL x 88.4 = μmol/L

**REFERENCE INTERVALS (NORMAL VALUES)  
(Based on rules CLSI C28-P Document)\***

**Serum / Plasma:**

Serum/plasma samples:

Men : 0.7 - 1.2 mg/dl (62 - 105 μmol/l)

Women: 0.6 - 1.1 mg/dl (53 - 97 μmol/l)

24h urine:

Men : 1000-2000 mg/24h (8.85-17.70 mmol/24h)

Women: 800-1800 mg/24h (7.08-15.93 mmol/24h)

\*It is recommended that each laboratory establish its own normal range.

**QUALITY CONTROL AND CALIBRATION**

All control sera with Creatinine values determined by this method can be used. We recommend:

Assayed Control Serum Normal

**A3910 ARCON N**

Assayed Control Serum Abnormal

**A3920 ARCON P**

This assay requires the use of a Creatinine Standard or a Creatinine Calibrator. We recommend:

**A2160S** Creatinine Standard (Conc. 2 mg/dl)

**A39050** Calibrator (ARCAL AUTO)

Any commercially available Standard or Calibrator suitable for this method may be used.

\*Calibration Stability: It is strongly depend of application to auto analyzers and auto analyzers specification. Calibration stability is 1 day in general.

\*Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

Quality control is recommended every morning. Calibration is not recommended if QC control values are acceptable. Reagent should be calibrated after lot changes.

**PERFORMANCE CHARACTERISTICS**

**Low Linearity:** 0.20 mg/dL

**High Linearity:** The method is linear up to 20 mg/dL.

If the value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Precision Studies (Based on CLSI EP5 Doc.):**

**Repeatability (within run)(intra-assay)**

mean (mg/dl)	SD (mg/dl)	CV%	n
1.25	0.03	2.60	10
3.87	0.07	1.90	10

**Reproducibility (run to run)(inter-assay):**

mean (mg/dL)	SD (mg/dL)	CV%	n
1.31	0.04	2.90	20
3.80	0.14	3.80	20

**Sensitivity (LOD) (Based on CLSI EP17 document):** The limit of detection is 0.2 mg/dl.

**Trueness:** A comparison between Archem and a commercially available product gave the following results:

Cretinine Archem = x

Creatinine competitor = y

n = 104

y = 0.982x - 0.081 mg/dl r<sup>2</sup> = 0.94

**Interference:** No interference was observed by the presence of:

Hemoglobin ≤ 200 mg/dL

Bilirubin ≤ 7 mg/dL

Lipids ≤ 600 mg/dL

These performance characteristics have been obtained using an analyzer. Results may vary if a different instrument or a manual procedure is used.

**NOTES**

1. For in vitro diagnostic use only. Do not pipette by mouth. Avoid contact with skin and mucous membranes.
2. All the calibrators, controls and some reagents must be considered as human&animal sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.
3. Material safety data sheet will be supplied on request.
4. Exercise the normal precautions required for handling laboratory reagents.
5. After measurements are taken, reagent bottles should capped and kept at 2-8°C. Caps of the reagents bottles can not be used between two different kind of reagent and between R1 & R2.
6. Reagents with different lot numbers should not be interchanged or mixed.
7. The linearity limit depends on the sample to reagent ratio.

**PRECAUTIONS AND WASTE DISPOSAL**

This product is made to be used in professional laboratories and by professional operators. Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.



R36/38 : Irritating to eyes and skin.  
 S20/21 : When using, do not eat, drink or smoke  
 S26 : In case of contact with eyes, rinse immediately with plenty of water and seek medical advice  
 S28 : After contact with skin wash immediately with plenty of water.  
 S36/37/39: Wear suitable protective clothing, gloves and eye/face protection.  
 S45 : In case of accident or if you feel unwell, seek medical advice immediately.  
 S56 : Dispose of this material and its container at hazardous or special waste collection point.  
 S57 : Use appropriate container to avoid environmental contamination.  
 S61 : Avoid release in environment. Refer to special instructions/safety data sheets  
 Please consult local regulations for a correct waste disposal.

### ABBREVIATIONS

CLSI : Clinical and Laboratory Standards Institute  
 CV% : Coefficient of Variation Percentage  
 EP : Evaluation Protocols  
 GLP : Good Laboratory Practice  
 IU : International Unit  
 mA : milliabsorbance  
 mL : milliliter  
 NCCLS: National Committee for Clinical Laboratory Standards  
 QC : Quality Control

### REFERENCES

1. Tietz Textbook of Clinical Chemistry, Second Edition, Burtis- Ashwood (1994).
2. Clinical and Laboratory Standards Institute (formerly NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. NCCLS Document EP05-A2.
3. HU Bergmeyer - Methods of enzymatic analysis, (1987).
4. Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 3rd ed. Washington: AACC Press; 1990.
5. Jaffe MZ. Ueber den Niederschlag, welchen Pikrinsäure in normalem Harn erzeugt and ueber eine Reaction des Kreatinins. *Zeitschrift Fuer Physiologische Chemie*, 1886;10:391-400 (Ger).
6. Hare, R.S., Endogenous Creatinine in Serum and Urine.
7. HENRY J.B., Clinical Diagnosis and Management; 17<sup>th</sup> edition, Saunder Publisher (1984)
8. Tietz NW. *Clinical Guide to Laboratory Tests*.

### SYMBOLS

<b>IVD</b>	Only for invitro diagnostic use
<b>LOT</b>	Lot of manufacturing
<b>R1</b>	Reagent 1
<b>R2</b>	Reagent 2
<b>CONC</b>	Concentration
<b>INGRED</b>	Reagent Ingredients
<b>SN</b>	Serial Number
<b>REF</b>	Reference Number (Catalog No)



Expiration date



Storage temperature interval



Read the directions



Biological risk



**Archem Diagnostics Industry LTD. ŞTİ.**  
 Organize Sanayi Bölgesi, Mutsan Sanayi Sitesi  
 M8 Blok No: 48 Başakşehir / İSTANBUL TURKEY  
 Tlf: + 90 212 444 08 92  
 Fax: +90 212 629 98 89  
[info@archem.com.tr](mailto:info@archem.com.tr)  
[www.archem.com.tr](http://www.archem.com.tr)