

CALCIUM

ARSENAZO III

Diagnostic reagent for determination of Calcium concentration.

Liquid. Mono Reagents. Store at 15°C - 25°C. For in Vitro Diagnostic Use. Do not freeze

Ref No	Pack	Ref No	Pack	Ref No	Pack	Ref No	Pack
A2060	5 x 100 mL	R2061	962 Tests	BY2060	6364 Tests	M2060	2273 Tests
A2061	5 x 50 mL	S2060	1920 Tests	BY2061	4545 Tests	M2061	909 Tests
A2062	5 x 25 mL	S2061	1440 Tests	BY2062	3182 Tests	L2060	5000 Tests
T2060	5294 Tests	K2060	3500 Tests	N2060	758 Tests	L2061	2667 Tests
T2061	2353 Tests	K2061	2500 Tests	N2061	424 Tests	DM2060	990 Tests
R2060	2885 Tests						

INTENDED USE

The test is applied for the quantitative determination of calcium in serum, plasma and urine.

TEST PRINCIPLE

Arsenazo (III) combines with calcium at slight acidic pH to form a blue complex, the absorbance of which is measured at 660 nm. The reaction has high specificity and interference from magnesium is avoided, due to pH.

For bichromatic analyzers, the reference wavelength must be set at 700 nm.

TEST PARAMETERS

Method	: Colorimetric, Endpoint, Increasing Reaction
Wavelength	: 660 nm (650 - 660 nm)
Temperature	: Room temperature, 37°C
Sample	: Serum or heparinized Plasma, acidified Urine
Linearity	: 1.5 mg/dL - 20 mg/dL (3.75 mmol/L)

REAGENT COMPOSITION

Arsenazo (III)	≤ 0.2 mM,
Good's buffer	≤ 50 mM
pH 6.8,	
Stabilizers.	

REAGENT PREPARATION

Reagents are ready to use.

REAGENT STABILITY AND STORAGE

Reagents are stable up to expiration date on labels at 15-25°C.

Once opened vials (reagent 1) are stable minimum 60 days at 15°C - 25°C at optimum conditions. On board stability is strongly related to auto analyzers cooling specification and carry-over values.

SAMPLE

Serum (preferred), plasma heparinate are collected by standard procedure. Do not use citrate, oxalate and EDTA as anticoagulant. Total calcium is stable 7 days at 2-8°C and for several months when frozen at -20°C.

Urine specimens should be collected in 20 to 30 ml of HCl 6M per 24/h specimen (1-2 ml for random urine) in order to prevent calcium salt precipitation.

Dilute sample urine 1:2 with redistilled water and multiply results by two.

TEST PROCEDURE

Sample Start

There have many ready application procedures dedicated to different kind of photometers and ready manual working process can be supplied on request.

There have many ready application procedures dedicated to different kind of biochemistry auto analyzers can be supplied on request.

CALCULATION

$$\frac{\text{Absorbance Sample}}{\text{Absorbance Standard}} \times \text{Conc. of Standard}$$

$$= \text{mg/dL (mmol/L) Calcium in sample}$$

Unit Conversion

$$\text{mg/dL} \times 0.2495 = \text{mmol/L}$$

**REFERENCE INTERVAL (NORMAL VALUES)
 (Based on CLSI C28-P Document)***

Serum/Plasma : 8.5 - 10.5 mg/dL
 (2.02 - 2.60 mmol/L)

Urine : 100-250 mg/24 Hour

*It is recommended that each laboratory establish its own reference range.

QUALITY CONTROL AND CALIBRATION

All control sera with calcium values determined by this method can be used. We recommend:

"ARCON N", Assayed Control Serum Normal
Cat.No. A3910

"ARCON P", Assayed Control Serum Abnormal
Cat.No. A3920

The assay requires the use of a Calcium Standard or a Calcium Calibrator. We recommend:

ARCHEM Standard
Cat.No. ASX2060S conc. 10 mg/dL

Any commercially available Standard or Calibrator suitable for this method may be used.

*Calibration Stability: It is strongly depend of application to auto analyzers and auto analyzers specification. Calibration stability is 15 days.

*Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

Quality control is recommended every morning. Calibration is not recommended if QC control values are acceptable. Reagent should be calibrated after lot changes.

PERFORMANCE CHARACTERISTICS

Low linearity: 1.5 mg/dL.

High Linearity: The assay is linear up to 20 mg/dL.

Above this concentration, dilute the sample with redistilled water or %0,9 NaCl and reassay multiplying the result by the dilution factor.

Linearity may considerably vary depending on the instrument used.

Precision Studies (Based on CLSI EP5 Doc.):
Repeatability (within run) (Intra-assay)

Mean conc.	SD	CV	n
8.89 mg/dL	0.10	1.1%	10
13.74 mg/dL	0.16	1.2%	10

Reproducibility (run to run) (Inter-assay)

Mean conc.	SD	CV	n
9.22 mg/dL	0.19	2.1%	10
14.04 mg/dL	0.23	1.7%	10

Sensitivity (LOD) (Based on CLSI EP17 document): The limit of detection is 0.1 mg/dL.

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

Interferences: No interference was observed by the presence of:

Hemoglobin ≤ 300 mg/dL
 Bilirubin ≤ 40 mg/dL

Lipids interferences are possible performing readings at single wavelength of 660 nm. To avoid interferences, perform a bichromatic reading at 660/700 nm.

Methods comparison: A comparison between Archem and a commercially available product gave the following results:

Calcium Archem = x
 Calcium competitor = y
 n = 97
 $y = 0.98x + 0.17 \text{ mg/dl} \quad r^2 = 0.94$

These performance characteristics have been obtained using an analyzer. Results may vary if a different instrument or a manual procedure is used.

NOTES

1. For in vitro diagnostic use only. Do not pipette by mouth. Avoid contact with skin and mucous membranes.
2. All the calibrators, controls and some reagents must be considered as human & animal sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.
3. Material safety data sheet will be supplied on request.
4. Exercise the normal precautions required for handling laboratory reagents.
5. After measurements are taken, reagent bottles should cap and kept at 2-8°C. Caps of the reagents bottles cannot be used between two different kind of reagent and between R1 & R2.
6. Reagents with different lot numbers should not be interchanged or mixed.
7. The linearity limit depends on the sample to reagent ratio.

PRECAUTIONS AND WASTE DISPOSAL

This product is made to be used in professional laboratories and by professional operators. Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

R36/38 : Irritating to eyes and skin.

S20/21 : When using, do not eat, drink or smoke.

S26 : In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S28 : After contact with skin wash immediately with plenty of water.

S36/37/39: Wear suitable protective clothing, gloves and eye/face protection.

S45 : In case of accident or if you feel unwell, seek medical advice immediately.

S56 : Dispose of this material and its container at hazardous or special waste collection point.

S57 : Use appropriate container to avoid environmental contamination.

S61 : Avoid release in environment. Refer to special instructions/safety data sheets.

Please consult local regulations for a correct waste disposal.

ABBREVIATIONS

CLSI : Clinical and Laboratory Standards Institute

CV% : Coefficient of Variation Percentage

EP : Evaluation Protocols

GLP : Good Laboratory Practice

IU : International Unit

mA : miliabsorbance

mL : milliliter

NCCLS: National Committee for Clinical Laboratory Standards

QC : Quality Control

REFERENCES

1. Young DS. Effects of Drugs on Clinical Laboratory Tests, AACC Press, Third Edition, Washington (1990).
2. Tietz NW. Clinical Guide to Laboratory Standards. 3rd ed. Philadelphia, PA: WB Saunders Company; 1995:102-105.
3. Tietz Textbook of Clinical Chemistry, Second Edition, Burtis- Ashwood (1994).
4. Endres DB, Rude RK. Mineral and bone metabolism. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1395-1457
5. Gitelman HJ. An improved procedure for the determination of calcium in biochemical specimens. Ana Biochem. 1967;18:521-531.
6. Clinical Chemistry vol. 38 N. 6 - 904-908 - (1992).
7. Clinical and Laboratory Standards Institute (formerly NCCLS). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. NCCLS Document EP05-A2.
8. Zak B., Epstein E., Babinski E.S., Review Calcium Methodologies, Annals of Clinical and Laboratory Science 5, 195-212 (1975).

SYMBOLS

IVD	Only for invitro diagnostic use
LOT	Lot of manufacturing
R1	Reagent 1
CONC	Concentration
INGRED	Reagent Ingredients
REF	Reference Number (Catalog No)
SN	Serial Number



Expiration date



Storage temperature interval



Read the directions



Biological risk



Archem Diagnostics Industry LTD. ŞTİ.

Organize Sanayi Bölgesi, Mutsan Sanayi Sitesi
 M8 Blok No: 48 Başakşehir / İSTANBUL TURKEY

Tlf: + 90 212 444 08 92

Fax: +90 212 629 98 89

info@archem.com.tr

www.archem.com.tr