

## CH CLEANER

Liquid. Store at +2/+30°C. For in Vitro Diagnostic Use (IVD). Do not freeze.

Produced for Siemens Atellica CH930.

Ref No	Package
CLNR10	2 x 1.5 L

### INTENDED USE

CH Cleaner is intended for in vitro diagnostic use in cleaning the dilution probe using the Atellica Chemistry Analyser.

### ACTIVE COMPOSITION

Sodium hypochlorite.....5%

### PREPARATION

CH Cleaner is ready for use.

### STABILITY AND STORAGE

Unopened CH Cleaner packages are stable at +2/+30 °C until the expiration date on the label. The onboard stability of opened CH Cleaner packages is 90 days. Dispose of products whose onboard stability has expired by removing them from the analyzer.

### TEST PRINCIPLE

CH Cleaner should be used as specified in the Atellica analyzers user manual.

### WARNINGS AND PRECAUTIONS

IVD: For in Vitro Diagnostic use only.

Do not use expired reagents.

Reagents with two different lot numbers should not be interchanged.

For professional use.

Follow Good Laboratory Practice (GLP) guidelines.



**H290:** May be corrosive to metals.

**H315:** Causes skin irritation.

**H318:** Causes serious eye damage.

**H400:** Very toxic to aquatic life.

**H411:** Toxic to aquatic life with long lasting effects.

**P234:** Keep only in original container.

**P264:** Cleaner hands thoroughly after handling.

**P280:** Wear protective gloves/protective clothing/eye protection/face protection.

**P273:** Avoid release to the environment.

**P305+P351+P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**P310:** Immediately call a POISON CENTER or doctor/physician.

**P390:** Absorb spillage to prevent material damage.

**P391:** Collect spillage.

**P501:** Dispose of contents and container in accordance with all local, regional, and national regulations.

### SYMBOLS

IVD

In Vitro Diagnostic Medical Device

LOT

Lot number

GTIN

Global Trade Item Number

REF

Reference Number

FOR USE WITH

Identifies products to be used together

PRODUCT OF TURKEY

Product of Turkey



Manufacturer



Expiration date



Temperature limitation



Humidity Limitation



Consult instructions for use



**Archem Sağlık Sanayi ve Tic. A.Ş.**

Mahmutbey Mah. Halkalı Cad. No:124 Kat:4

Bağcılar/İstanbul/Turkey

Tel: + 90 212 444 08 92

Fax: +90 212 629 98 89

info@archem.com.tr www.archem.com.tr

