

CRP WR Turbidimetric Calibrator

REF 01R91-01

FOR USE WITH
ARCHITECT

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package.

INTENDED USE

CRP WR Turbidimetric Calibrator is for use in the calibration of the CRP WR Turbidimetric assay on the Architect cSystems.

CRP WR Turbidimetric Calibrator

REF: 01R91-01

Package: 5 x 1.5 mL **Liquid**

For use with;

CRP WR Turbidimetric reagent

REF: *02R04-31

REF: *02R04-21

CONTENTS / MATERIALS PROVIDED

CRP WR Turbidimetric Calibrator (REF: 01R91-01) is a bovine serum based solution.

Sodium azide (0.09 %) is added as a preservative.

MATERIALS REQUIRED BUT NOT PROVIDED

1. CRP WR Turbidimetric reagent
REF: *02R04-31
REF: *02R04-21

STANDARDIZATION

The traceability of the method is verified using European Reference Materials (ERM) DA 472/IFCC or 85/506 NIBSC for CRP WR Calibrator.

PRECAUTIONS (General and Users)

IVD

For *in Vitro* Diagnostic Use.

Do not use components beyond the expiration date.

Contains sodium azide.

EUH032 Contact with acids liberates very toxic gas.

Do not mix materials from different kit lot numbers.

Animal/Plasma Cell sourced material.

Calibrator can contain preservatives (as Sodium Azide or others) which total concentration is lower than the limits mentioned in Directive 67/548/ CEE and 88/379 CEE.



Treat as potentially infectious.

Safety Data Sheets are available at www.archem.com.tr or contact your local representative. Emergency Phone Number: In emergency situation please contact the National Poison

Solidarity Center (Turkey) for information in Turkish:(114) (24 hours/7 days). Please contact with Archem for information in English: +(90) 212 444 0 892 (Turkey) (8:30-18:30/ except Sundays)

STORAGE AND STABILITY

1. Unopened CRP WR Turbidimetric Calibrator is stable until the expiration date when stored +2 to +8 °C.
2. Opened calibrators are stable **45 days** when stored tightly capped in a dark place +2 to +8°C and kept free of contamination.

PREPARATION OF CALIBRATOR

Liquid, ready to use. No need for preparation.

INDICATIONS OF INSTABILITY OR DETERIORATION

Presence of extreme turbidity or microbial growth may indicate deterioration.

REFERENCES

1. Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 4th ed. St. Louis , MO, Elsevier Saunders; 2006:2263.
2. Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR Part 1910.1030). Fed. Register. July 1, 2001; 17:260–273.
3. Directive 2000/54/EC. Official Journal of the European Communities No. L262 from October 17, 2000.
4. Tietz, N. W. (1983) *Clinical guide to Laboratory Tests*, p384 W.B. Saunders Co., Philadelphia.







TRADEMARKS

CRP WR Turbidimetric Calibrator is a trademark of ARCHEM Diagnostic A.Ş. In various jurisdictions. The ARCHITECT c System family of instruments consists of c 4000, c8000 and c 16000 instruments. ARCHITECT, c4000, c8000, c16000, c System are trademark of Abbott Laboratories in various jurisdictions All trademarks are property of their respective owner(s).

Customer Service:

Contact Archem Diagnostics Industry A.Ş.

**CRP WR Turbidimetric
Calibrator****SYMBOLS**

IVD	In Vitro Diagnostic Medical Device
LOT	Lot number
CAL	Calibrator
GTIN	Global Trade Item Number
REF	Reference Number
FOR USE WITH	Identifies products to be used together
PRODUCT OF TURKEY	Product of Turkey
	Manufacturer
	Expiration date
	Temperature limitation
	Consult instructions for use
	Caution
	Sufficient for

chem
GNOSTICS

Archem Sağlık Sanayi ve Tic. A.Ş.
Mahmutbey Mah. Halkalı Cad. No:124 Kat:4
Bağcılar/İstanbul/Turkey
Tel: + 90 212 444 08 92
Fax: +90 212 629 98 89
info@archem.com.tr www.archem.com.tr

