TÜRKÇE (TR)

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ELÜENT 80A, 80B

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For use with ADAMS™ A1C

ELUENT 80A, 80B

Always keep this package insert on hand and carefully read it before use.

INTENDED USE:

These products are intended for the determination of Hemoglobin A1C (HbA1C) in human whole blood for use with ADAMS™ A1C HA-8180 series. For in vitro diagnostics use only.

SUMMARY AND PRINCIPLE:

The measurement of stable HbA1C is very useful in the diagnosis and treatment of diabetes and observations of diabetic conditions. The HA-8180 measures HbA1C in blood using reversed-phase cation exchange chromatography. Blood sample into several hemoglobin components based on high performance liquid chromatography (HPLC).

WARNINGS AND PRECAUTIONS

Warnings:

- As each eluent contains an oxidizer, if eluent spills onto the countertop of floor, immediately and carefully wipe it up with a cloth, then rinse the collected liquid from the cloth with plenty of water. Then, rinse the collected liquid from the cloth with plenty of water.
- If the spilled liquid dries and crystallizes, wipe up the spill using a cloth moistened with water. Then, rinse the collected liquid from the cloth with plenty of water.
- Never contact spilled liquid or crystals with reducing agents such as alcohol or ascorbic acid.
- Dilute any remaining reagent in eluent packs with plenty of water before discarding.
- Each eluent contains less than 0.06% mucosal-irritative sodium azide (NaN₃) as a preservative.

- Sodium azide may also react with copper or lead whereby potentially generating explosive metal azide.
- Sodium azide may also react with acid leading to the generation of toxic explosive hydrogen azide.
- Therefore, when these materials are disposed of, you should flush them with an adequate volume of water in order to prevent any residual reagent in the water service pipe.
- If any of these eluents make contact with eyes or mouth, immediately wash with plenty of water and consult a doctor.
- If it makes contact with skin, wash with plenty of water.
- Discard liquid waste than contains blood sample in accordance with local regulations for biohazardous waste.

Precautions:

- For in vitro diagnostic use only.
- These products are only for ADAMS™ A1C HA-8180 series. DO NOT attempt to use with other instruments.
- If these products are stored in a cold place, such as a refrigerator, allow them to adjust to the same environment as the instrument for at least one hour before placing them on the instrument.
- DO NOT use eluent packs beyond their expiration dates.
- The expiration dates are written on both the box and pack label.
- Place each eluent pack in the specified position on the bottle tray.
- Attach the proper bottle cap with nozzle to each Eluent A pack and Eluent B pack.
- Connecting the wrong bottle cap to the wrong pack may cause mixing of eluents, producing inaccurate measurement results.
- Replace the eluent pack with a new one even if a small volume of eluent remains.
- Adding new eluent to the old pack may cause inaccurate measurement results.
- Always handle these eluents with care.
- Avoid physical shock.

MEASUREMENT PROCEDURE

Refer to the operating manual for ADAMS™ A1C HA-8180 series.

STORAGE:

- Store unopened eluent packs at a temperature between 3°C and 30°C, avoiding direct sunlight.
- Once opened, use the eluent within one month, even if it is before expiration date.
- If the instrument is not to be used for one week or more, remove the eluent packs from the instrument. Then, cap them tightly and store them at a temperature between 3°C and 30°C in order to prevent liquid degradation and/or concentration.

MATERIALS PROVIDED:



1. ELUENT 80A: 600 mL X 4 packs/ box, phosphate buffer
2. ELUENT 80B: 600 mL X 2 packs/ box, phosphate buffer

Materials of the pack:

Aluminum pack with PP (polypropylene) cap
The aluminum pack consists of polyethylene (PE), aluminum foil and polyethylene terephthalate (PET).

 Archem Sağlık Sanayi ve Tic. A.Ş.
Mahmutbey Mah. Halkalı Cad. No:124 Kat:4
Bağcılar/İstanbul/Türkiye
Tel: +90 212 444 08 92
Fax: +90 212 629 98 89
info@archem.com.tr www.archem.com.tr

SYMBOLS	
	In Vitro Diagnostic Medical Device
	Lot Number
	Global Trade Item Number
	Reference Number
	Identifies Products to Be Used Together
	Product of Turkey
	Manufacturer
	Expiration Date
	Temperature Limits
	Consult Instructions for Use