TÜRKÇE (TR)

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HEMOLYSIS WASHING SOLUTION 80H

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For use with ADAMS™ A1C

HEMOLYSIS WASHING SOLUTION 80H

Always keep this package insert on hand and carefully read it before use.

INTENDED USE:

These products are intended for the determination of hemoglobin A1c (HbA1c) in human whole blood for use with ADAMS™ A1C HA-8180 series. Use this product to dilute whole blood samples and to wash the tubes in which sample flows. For in vitro diagnostics use only.

SUMMARY AND PRINCIPLE:

The measurement of stable HbA1c is very useful in the diagnosis and treatment of diabetes and observations of diabetic conditions. The HA-8180 measures HbA1c in blood using reversed-phase cation exchange chromatography. Blood sample diluted with the hemolysis washing solution is sent to the column, which fractionates the sample into several hemoglobin components based on high performance liquid chromatography (HPLC).

WARNINGS AND PRECAUTIONS

Warnings:

- This product contains less than 0.02% mucosal-irritative sodium azide (NaN₃) as a preservative.
- Sodium azide may also react with copper or lead whereby potentially generating explosive metal azide.
- Sodium azide may also react with acid leading to the generation of toxic explosive hydrogen azide.
- Therefore, when these materials are disposed of, you should flush them with an adequate volume of water in order to prevent any residual reagent in the

water servicepipe.

- If this product makes contact with eyes or mouth, immediately wash with plenty of water and consult a doctor.
- If it makes contact with skin, wash with plenty of water
- Discard liquid waste than contains blood sample in accordance with local regulations for biohazardous waste.

Precautions:

- For in vitro diagnostic use only.
 - This product is only for ADAMS™ A1C HA-8180 series. DO NOT attempt to use with other instruments.
 - If this product is stored in a cold place, such as a refrigerator, allow it to adjust to the same environment as the instrument for at least one hour before placing it on the instrument.
 - DO NOT use hemolysis washing solution bottle beyond its expiration dates.
 - The expiration date is written on both the box and bottle label.
 - Place the hemolysis washing solution bottle in the specified position on the bottle tray and correctly attach the bottle cap with the nozzle for the hemolysis washing solution.
 - Replace the hemolysis washing solution bottle with a new one even if a small volume of solution remains.
 - Adding new solution to the old bottle may cause inaccurate measurement results.
 - Always handle this product with care.
 - Avoid physical shock.
- MEASUREMENT PROCEDURE**
Refer to the operating manual for ADAMS™ A1C HA-8180 series.

STORAGE:

- Store unopened hemolysis washing solution bottles at a temperature between 3°C and 30°C, avoiding direct sunlight.
- Once opened, use the hemolysis washing solution within one month, even if it is before expiration date.
- If the instrument is not to be used for used for one week or more, remove the hemolysis washing solution bottle from the instrument. Then, cap them tightly and store them at a temperature between 3°C and 30°C in order to prevent liquid degradation and/or concentration.

MATERIALS PROVIDED:

HEMOLYSIS WASHING SOLUTION 80H: 2 L X 3 packs/ box, phosphate buffer containing surfactant.




Materials of the pack:

Polyethylene (PE) bottle with polypropylene (PP) cap.

 Archem Sağlık Sanayi ve Tic. A.Ş.
Mahmutbey Mah. Halkalı Cad. No:124 Kat:4
Bağcılar/İstanbul/Türkiye
Tel: +90 212 444 08 92
Fax: +90 212 629 98 89
info@archem.com.tr www.archem.com.tr

SYMBOLS

	In Vitro Diagnostic Medical Device
	Lot Number
	Global Trade Item Number
	Reference Number
	Identifies Products to Be Used Together
	Product of Turkey
	Manufacturer
	Expiration Date
	Temperature Limits
	Consult Instructions for Use