

CONCENTRATED CLEANER

Ref No	Packing	Ref No	Packing
5141M	50 mL	5940MN	500 mL

Changes made in the instructions for use are marked as grey

INTENDED USE

It is used to clean the instrument regularly.

TEST PRINCIPLE

For detailed usage instruction, please refer to the Operator's Manual of the instruments.

REAGENT COMPONENTS

Surfactant	≤0.2%
Sodium hypochlorite	≤12%
Sodium hydroxide	≤5%

REAGENT PREPARATION

Reagents are ready for use.

REAGENT STORAGE AND STABILITY

- It is important to store in dark place and keep container closed when not in use.
- The product will be stable for 18 months when stored at 2°C to 30°C (35.6 °F to 86 °F) in a dark place, with relative humidity not above 90%.
- The working temperature range of the product is consistent with that of its applicable instruments.
- The open-vial validity is 60 days.

SAMPLE REQUIREMENTS

Not applicable.

REFERENCE INTERVALS / MEDICAL DECISION LEVELS

Not applicable.

PERFORMANCE CHARACTERISTICS

Not applicable.

RESULT ELABORATION

Not applicable.

LIMITATIONS

Not applicable.

WARNINGS AND PRECAUTIONS

IVD: For in Vitro Diagnostic use only.
 Do not use expired reagents.
 Reagents with two different lot numbers should not be interchanged.
 For professional use.
 Follow Good Laboratory Practice (GLP) guidelines.

CAUTION: Human source samples are processed with this product. All human source samples must be treated as potentially infectious materials and must be handled in accordance with OSHA standards.

Danger

H411	Toxic to aquatic life with long lasting effects.
H400	Very toxic to aquatic life.
H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
EUH031	Contact with acids liberates toxic gas.
EUH206	Warning!

Do not use together with other products. May release dangerous gases (chlorine).

Prevention

P260	Do not breathe mist/vapours/spray.
P264	Wash all exposed external body areas thoroughly after handling.
P280	Wear protective gloves, protective clothing, eye protection and face protection.
P234	Keep only in original packaging.
P273	Avoid release to the environment.

Response

P301+P330+P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER/ doctor /physician/first aider.
P363	Wash contaminated clothing before reuse.
P390	Absorb spillage to prevent material damage.

P391 Collect spillage.

P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Storage

P405 Store locked up.

Disposal



P501 Dispose the vials and contents according to the local regulations.

REFERENCES







Not applicable.



Archem Sağlık Sanayi ve Tic. A.Ş.
Mahmutbey Mah. Halkalı Cad. No:124 Kat:4
Bağcılar/Istanbul/Türkiye
Tel: + 90 212 444 08 92
Fax: +90 212 629 98 89
info@archem.com.tr www.archem.com.tr



SYMBOLS

IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
R1	Reagent 1
R2	Reagent 2
GTIN	Global Trade Item Number
REF	Reference Number
GLP	Good Laboratory Practices
FOR USE WITH	Identifies Products to Be Used Together
PRODUCT OF TURKEY	Product of Turkey
	Manufacturer
	Expiration Date
	Temperature Limits
	Consult Instructions for Use
	Caution
	Number of Tests