

CUVETTE WASH

Liquid. Store at +5/+25°C. For in Vitro Diagnostic Use (IVD). Do not freeze.
Produced for Siemens Advia Chemistry Analyser.

Ref No	Package
AD3000	2L

Changes made in the instructions for use are marked as grey.

INTENDED USE

Cuvette Wash is intended for in vitro diagnostic use in washing cuvettes using the Advia Chemistry Analyser.

ACTIVE COMPOSITION

Sodium hydroxide.....4%

PREPARATION

Cuvette Wash is ready for use.

STABILITY AND STORAGE

Unopened Cuvette Wash packages are stable at +5/+25 °C until the expiration date on the label. The onboard stability of opened Cuvette Wash packages is 90 days. Dispose of products whose onboard stability has expired by removing them from the analyzer.

TEST PRINCIPLE

Cuvette Wash should be used as specified in the Advia Chemistry Analyser user manual.





WARNINGS AND PRECAUTIONS

IVD: For in Vitro Diagnostic use only.
Do not use expired reagents.
Reagents with two different lot numbers should not be interchanged.
For professional use.
Follow Good Laboratory Practice (GLP) guidelines.



H290: May be corrosive to metals.
H314: Causes severe skin burns and eye damage.
P234: Keep only in original container.
P264: Wash hands thoroughly after handling.
P280: Wear protective gloves/protective clothing/eye protection/face protection.
P301+P330+P331: IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
P303+P361+P353: IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
P310: Immediately call a POISON CENTER or doctor/physician.

P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P390: Absorb spillage to prevent material damage.
P501: Dispose of contents and container in accordance with all local, regional, and national regulations.

SYMBOLS	
IVD	In Vitro Diagnostic Medical Device
LOT	Lot number
GTIN	Global Trade Item Number
REF	Reference Number
FOR USE WITH	Identifies products to be used together
PRODUCT OF TURKEY	Product of Turkey
	Manufacturer
	Expiration date
	Temperature limitation
	Consult instructions for use



Archem Sağlık Sanayi ve Tic. A.Ş.
Mahmutbey Mah. Halkalı Cad. No:124 Kat:4
Bağcılar/İstanbul/Turkey
Tel: + 90 212 444 08 92
Fax: +90 212 629 98 89
info@archem.com.tr www.archem.com.tr

